



Total Alkalinity Test Strips, 0-240 mg/L, 50 tests

Termékszám:	2744850
HUF Ár (nettó):	Elérhetőségünk
Elérhető	

Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to drinking water, food, and beverage processors; boiler operators; and aquatic biologists. With Hach's alkalinity test strips, technicians in the field can test many samples in only a few minutes, and make immediate evaluations on-site, without any measuring, set-up, clean-up, or chemical handling. Our alkalinity test strips are also used in laboratories all over the world for pre-test screening-to detect the presence of materials that might interfere with lab testing. Contains 50 tests.

Easy to use, disposable, and inexpensive

A great way to obtain quick, quantitative answers in the field or in the lab

Műszaki adatok

Egység:	50 tests
Legkisebb beosztás:	Steps: 0, 40, 80, 120, 180, 240 mg/L (ppm)
Méréstartomány:	0 - 240 mg/L
Paraméter:	Total Alkalinity - Phenolphthalein
Platform:	Tesz csík
Tesztek száma:	50
Tesztkészlet model:	Test Strips