



## Digital Colony Counter Pen

<b>Termékszám:</b>	<b>25045000</b>
HUF Ár (nettó):	Elérhetőségünk
Szállítás 2 héten belül	

Counter Pen is a combination marker and digital counter that can be used to count cell growth on petri dishes, inventory parts, and tally samples received. To ensure an accurate tally, the counter beeps each time a mark is made, and the easy-to-read LCD displays the total (up to 99,999). Pen features on/off switch and comes with cap.

**Count cell growth on petri dishes**

**Inventory parts**

**Tally samples received**

---

### Műszaki adatok

Anyag:	Durable ABS plastic
Térfogat/ csomagméret:	1