



Be Right™



Beaker, polypropylene, low form, 50 mL, 12/pk

Termékszám:

108071

HUF Ár (nettó):

Elérhetőségünk

Szállítási dátum nem látható

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes $\geq 250\text{mL}$. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Műszaki adatok

Anyag:	Polypropylene, Low Form
Anyaga:	Polypropylene
Átmérő, külső:	51.1 mm with spout
Egység:	12 /pk
Kapacitás:	50 mL