



Be Right™



## Beaker, polypropylene, low form, 250 mL

<b>Termékszám:</b>	<b>108046</b>
HUF Ár (nettó):	Elérhetőségünk
Elérhető	

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes  $\geq 250$  mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

---

### Műszaki adatok

Anyag:	Polypropylene, Low Form
Anyaga:	Polypropylene
Átmérő, külső:	83.6 mm with spout
Egység:	1 each
Kapacitás:	250 mL