



**Be Right™**



## **Funnel, analytical, 65 mL approx. volume**

<b>Termékszám:</b>	<b>108367</b>
HUF Ár (nettó):	Elérhetőségünk
Szállítási dátum nem látható	

Standard long stem funnel meets the needs of a wide variety of applications. Mounded ribbing prevents "air lock" commonly occurring during funneling and filtrations. A steep, 58-degree cone slope makes for fast and easy flow of either liquids or small solids. Polypropylene construction is resistant to chemicals and makes this funnel fully autoclavable for sterilization.

---

### **Műszaki adatok**

Térfogat/ csomagméret: 1 each