



## Butterfield's buffered phosphate

<b>Termékszám:</b>	<b>2319110</b>
HUF Ár (nettó):	Elérhetőségünk
Szállítási dátum nem látható	

Butterfields Buffered Phosphate, 450 mL bottle. pk/10.

### **Complete Solutions for Testing Recreational Water**

Hach offers a variety of media and apparatus to meet all of your recreational testing needs.

### **Complete Solutions for Testing Wastewater**

Hach has prepared media and convenient apparatus to make wastewater and non-potable water testing fast and easy.

### **Complete Solutions for Testing Drinking Water**

Whichever method a customer chooses for examining drinking water samples, Hach has prepared media and convenient apparatus to make testing easy and fast.

---

## **Műszaki adatok**

Egység:	10/pk
Térfogat:	450 mL
Térfogat/ csomagméret:	10/pk