



Bottle rack, plastic-coated, 20 compartments

Termékszám:	2551801
HUF Ár (nettó):	Elérhetőségünk
Szállítási dátum nem látható	

Plastic-coated wire rack provides corrosion resistant transport system for sample bottles and bags. Racks feature 15 or 20 compartments with double-grid bottom to prevent loss of samples during transport. 15 compartment rack is ideal for Whirl-Pak® sample bags (4 oz). 20 compartment rack fits most bottles associated with standard drinking water sampling. Autoclavable.

Műszaki adatok

Speciális tulajdonság:	Rack, 20 Bottle Compartments
Térfogat/ csomagméret:	1