



Beaker, griffin, low form, glass, 50 mL, 12/pk

Termékszám:	50071H
HUF Ár (nettó):	Elérhetőségünk
Elérhető	

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologic

Műszaki adatok

Anyag:	Glass, Griffin, Low Form
Anyaga:	Glass
Egység:	12 /pk
Kapacitás:	50 mL